

Aortic catheters





BODY LENGTH: 21 CM


Description	Size	Standard	With vented red cap
With connector	20 Fr	AC201341	AC201341V
	22 Fr	AC221341	AC221341V
	24 Fr	AC241341	AC241341V
	26 Fr	AC261341	AC261341V
	30 Fr	AC301341	AC301341V
With luer connector	22 Fr	AC221342	AC221342V
	24 Fr	AC241342	AC242342V
	26 Fr	AC261342	AC261342V
	30 Fr	AC301342	AC301342V

PRODUCT INFORMATION

Andocor aortic catheters are intended for perfusion of the ascending aorta during cardiopulmonary bypass.

-  Connection site: 3/8"
Qty/box: 20 Pcs
-  Easy tip orientation thanks to indexing marker
Smooth bevelled tip
Moveable suture ring

AORTIC CATHETER PRESSURE DROP VS FLOW



Special configurations

HARD TIP

BODY LENGTH: 21 CM

Description	Size	Standard	With vented red cap
With connector	24 Fr	AC241341H	AC241341VH
With luer connector	24 Fr	AC241342H	AC241342VH

WITH SMALL RING, HARD TIP



BODY LENGTH: 21 CM

Description	Size	Standard	With vented red cap
With connector	24 Fr	-	AC241341VHS
With luer connector	24 Fr	-	AC241342VHS

HARD TIP, LONG



BODY LENGTH: 36 CM

Description	Size	With vent plug	With vented red cap
With luer connector	22 Fr	AC221342HL	AC221342VHL
	24 Fr	AC241342HL	AC241342VHL
With luer connector and side holes	22 Fr	AC221342HLW	AC221342VHLW
	24 Fr	AC241342HLW	AC241342VHLW