


Aortic root cannulae



Description	Tip size	Needle size	Reference
Type A Without vent line	6 Fr	16 G	ARC07
	7 Fr	14 G	ARC04
	9 Fr	12 G	ARC01
Type B Without vent line With suture ring	6 Fr	16 G	ARC07R
	7 Fr	14 G	ARC04R
	9 Fr	12 G	ARC01R
Type C With vent line	6 Fr	16 G	ARC08
	7 Fr	14 G	ARC05
	9 Fr	12 G	ARC02
Type D With vent line	6 Fr	16 G	ARC08S
	7 Fr	14 G	ARC05S
	9 Fr	12 G	ARC02S


PRODUCT INFORMATION

Andocor aortic root cannulae are intended for use during cardiopulmonary bypass for antegrade delivery of cardioplegia solutions and venting of the heart. The cannula may also be used to aspirate air from the aorta at the conclusion of the cardiac procedure.

Q Length: 15 cm
 Connection size vent line: 1/4"
 Qty/box: 20 Pcs

+ Easy of use
 Anti-bleeding valve helps reduce blood loss

AORTIC ROOT CANNULAE PRESSURE DROP VS FLOW




WITH 2 CLAMPS




Description	Tip size	Needle size	Reference
With vent line	7 Fr	14 G	ARC10

WITH ANTI-BEELDING VALVE



Description	Tip size	Needle size	Reference
Without vent line	6 Fr	16 G	ARC07V
	7 Fr	14 G	ARC04V
	9 Fr	12 G	ARC01V
With vent line	6 Fr	16 G	ARC08V
	7 Fr	14 G	ARC05V
	9 Fr	12 G	ARC02V

MINIMAL INVASIVE TYPE, LENGTH: 28 CM



Description	Tip size	Needle size	Standard	With anti-bleeding valve
Without vent line	7 Fr	14 G	ARC14	ARC14V
With vent line	7 Fr	14 G	ARC15S	-