

Arterial cannulae

REINFORCED, CURVED TIP WITH FLANGE



TYPE A - BODY LENGTH: 30 CM

TYPE B - BODY LENGTH: 21 CM

Description	Size	TYPE A - BODY LENGTH: 30 CM			TYPE B - BODY LENGTH: 21 CM		
		Standard	With vent plug	With vented red cap	Standard	With vent plug	With vented red cap
Without connector	20 Fr	A2016	A2016V	-	B2016	B2016V	-
	22 Fr	A2216	A2216V	-	B2216	B2216V	-
	24 Fr	A2416	A2416V	-	B2416	B2416V	-
With connector	20 Fr	A20161	-	A20161V	B20161	-	B20161V
	22 Fr	A22161	-	A22161V	B22161	-	B22161V
	24 Fr	A24161	-	A24161V	B24161	-	B24161V
With luer connector	20 Fr	A20162	-	A20162V	B20162	-	B20162V
	22 Fr	A22162	-	A22162V	B22162	-	B22162V
	24 Fr	A24162	-	A24162V	B24162	-	B24162V

PRODUCT INFORMATION


Andocor arterial cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.

All versions with a vented red cap can be equipped with a vent plug. Upon request the luer connector can be mounted at an angle 180° - facing up.

Q Connection site: 3/8"
Qty/box: 20 Pcs

+ Kink resistant reinforced tubing
Easy tip orientation thanks to indexing marker
Smooth plastic tip

ARTERIAL CANNULA PRESSURE DROP VS FLOW



Arterial cannulae

REINFORCED, STRAIGHT TIP WITH RING



TYPE A



TYPE B

TYPE A - BODY LENGTH: 30 CM

TYPE B - BODY LENGTH: 21 CM

Description	Size	TYPE A - BODY LENGTH: 30 CM			TYPE B - BODY LENGTH: 21 CM		
		Standard	With vent plug	With vented red cap	Standard	With vent plug	With vented red cap
Without connector	20 Fr	A2013	A2013V	-	B2013	B2013V	-
	22 Fr	A2213	A2213V	-	B2213	B2213V	-
	24 Fr	A2413	A2413V	-	B2413	B2413V	-
With connector	20 Fr	A20131	-	A20131V	B20131	-	B20131V
	22 Fr	A22131	-	A22131V	B22131	-	B22131V
	24 Fr	A24131	-	A24131V	B24131	-	B24131V
With luer connector	20 Fr	A20132	-	A20132V	B20132	-	B20132V
	22 Fr	A22132	-	A22132V	B22132	-	B22132V
	24 Fr	A24132	-	A24132V	B24132	-	B24132V

PRODUCT INFORMATION


Andocor arterial cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.

All versions with a vented red cap can be equipped with a vent plug. Upon request the luer connector can be mounted at an angle 180° - facing up.

Q Connection site: 3/8"
Qty/box: 20 Pcs

+ Kink resistant reinforced tubing
Easy tip orientation thanks to indexing marker
Smooth plastic tip

ARTERIAL CANNULA PRESSURE DROP VS FLOW



Arterial cannulae

NON-REINFORCED, CURVED TIP WITH FLANGE



BODY LENGTH: 22 CM

Description	Size	Standard	With vented red cap
Without connector	20 Fr	A2006	-
	22 Fr	A2206	-
	24 Fr	A2406	-
With connector	20 Fr	A20061	A20061V
	22 Fr	A22061	A22061V
	24 Fr	A24061	A24061V
With luer connector	20 Fr	A20062	A20062V
	22 Fr	A22062	A22062V
	24 Fr	A24062	A24062V

PRODUCT INFORMATION


Andocor arterial cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.

All versions with a vented red cap can be equipped with a vent plug. Upon request the luer connector can be mounted at an angle 180° - facing up.

Q Connection site: 3/8"
Qty/box: 20 Pcs


+ Kink resistant non-reinforced tubing
Easy tip orientation thanks to indexing marker
Smooth plastic tip

ARTERIAL CANNULA PRESSURE DROP VS FLOW



Arterial cannulae

NON-REINFORCED, STRAIGHT TIP WITH RING



BODY LENGTH: 22 CM

Description	Size	Standard	With vented red cap
Without connector	20 Fr	A2003	-
	22 Fr	A2203	-
	24 Fr	A2403	-
With connector	20 Fr	A20031	A20031V
	22 Fr	A22031	A22031V
	24 Fr	A24031	A24031V
With luer connector	20 Fr	A20032	A20032V
	22 Fr	A22032	A22032V
	24 Fr	A24032	A24032V

PRODUCT INFORMATION


Andocor arterial cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.

All versions with a vented red cap can be equipped with a vent plug. Upon request the luer connector can be mounted at an angle 180° - facing up.

Q Connection site: 3/8"
Qty/box: 20 Pcs

+ Kink resistant non-reinforced tubing
Easy tip orientation thanks to indexing marker
Smooth plastic tip

ARTERIAL CANNULA PRESSURE DROP VS FLOW



Arterial cannulae

NON-REINFORCED, ONE PIECE BODY



BODY LENGTH: 22 CM

Description	Size	Reference
Without connector, with suture ring	18 Fr	A1826
	20 Fr	A2026
	21 Fr	A2126
	22 Fr	A2226
	24 Fr	A2426

PRODUCT INFORMATION

The cannula consists of kink-resistant PVC tubing, it is constructed of a conical, single piece body terminating in a bevelled, curved tip. Each cannula is supplied with a fixed depth ring and has a blue indexing marker to allow precise positioning of the cannula.



Connection site: 3/8" - 18 Fr: 1/4"
Qty/box: 20 Pcs



Kink resistant pvc tubing
Easy tip orientation thanks to indexing marker

