


# Cardioplegia adapters



## STRAIGHT ADAPTERS

Description	Tubing length	Reference
1/4" / male luer connection	10 cm	CSY15
Female / male luer connection	10 cm	CSY16




## PERFUSION ADAPTERS

Description	Length	Reference
1/4" to male luer connector	7 cm	CSY17
3/8" to male luer connector	7 cm	CSY18




Description	Length	Reference
1/4" to male luer connector with stopcock	10 cm	CSY19



## PERFUSION Y-ADAPTERS

Description	Tubing length	Reference
Female luer lock in combination with 2 male luer lock connectors.	30 cm	CSY03
	35 cm	CSY03A

Description	Tubing length	Reference
Male luer lock in combination with 2 female luer lock connectors.	45 cm	CSY09
	30 cm	CSY20



PERFUSION SET WITH VENT LINE

Description	Tubing length	Reference
Male luer lock in combination with 1 female luer lock and 1 cone connector.	30 cm	CSY08
	18 cm	CSY13



PERFUSION RECIRCULATION ADAPTERS

Description	Tubing length	Reference
Male luer lock in combination with 2 cone connectors.	30 cm	CSY06
	13,5 cm	CSY25


Description	Tubing length	Reference
Female luer lock in combination with 3 male luer lock connectors	30 cm	CSY21

PRODUCT INFORMATION

Andocor cardioplegia adapters connect to the aortic root cannula or vessel cannula and are intended for the delivery of cardioplegia solution or venting the heart during cardiopulmonary bypass surgery.

 Qty/box: 20 Pcs


# Cardioplegia adapters



CSY04


## MULTIPLE PERFUSION ADAPTER FOR AORTIC ROOT AND THREE VEIN GRAFTS

Tubing length: 32 cm	Reference
Female luer lock in combination with 4 male luer lock connectors	CSY04




CSY05

Tubing length: 32 cm	Reference
Female luer lock in combination with 4 male luer lock connectors with 3 attached vessel cannula	CSY05



CSY24

Tubing length: 32 cm	Reference
Female luer lock in combination of 4 male luer lock connectors	CSY24



CSY23


Tubing length: 32 cm	Reference
Female luer lock in combination with 6 male luer lock connectors	CSY23



CSY12

## MULTIPLE PERFUSION ADAPTER FOR RETROGRADE AND AORTIC ROOT WITH VENT LINE


Tubing length: 45 cm	Reference
Female luer lock in combination with 2 male luer lock connectors and cone connector	CSY12



MULTIPLE PERFUSION ADAPTER WITH VENT LINE

Tubing length: 32 cm	Reference
Female luer lock in combination with 4 male luer lock connectors and cone connector	CSY22

Tubing length: 32 cm	Reference
Female luer lock in combination with 4 male luer lock connectors with 3 attached vessel cannula and cone connector	CSY07



MULTIPLE PERFUSION SET FOR THREE VEIN GRAFTS

Tubing length: 30 cm	Reference
Female luer lock in combination with 4 male luer lock connectors and cone connector	CSY10

Tubing length: 30 cm	Reference
Male luer lock in combination with 3 male luer lock connectors with attached vessel cannula	CSY11

PRODUCT INFORMATION

Andocor cardioplegia adapters connect to the aortic root cannula or vessel cannula and are intended for the delivery of cardioplegia solution or venting the heart during cardiopulmonary bypass surgery.

 Qty/box: 20 Pcs