

Flex line venous cannulae

REINFORCED WITH LIGHTHOUSE TIP



PRODUCT INFORMATION


Andocor venous cannulae are intended to drain venous blood from the right atrium and caval veins.

Q Connection site: 3/8" - 36 Fr: 1/2"
Length: 39 cm
Qty/box: 20 Pcs


+ Single piece construction
Open and flexible lighthouse tip
Easy positioning
Insertion depth marker

Size	Straight tip	Right angled 90° tip
20 Fr	01V201L8	01V201L9
22 Fr	01V221L8	01V221L9
24 Fr	01V241L8	01V241L9
26 Fr	01V261L8	01V261L9
28 Fr	01V281L8	01V281L9
30 Fr	01V301L8	01V301L9
32 Fr	01V321L8	01V321L9
34 Fr	01V341L8	01V341L9
36 Fr	01V361L8	01V361L9

VENOUS CANNULA PRESSURE DROP VS FLOW



VENOUS CANNULA PRESSURE DROP VS FLOW



Special configurations

CONNECTION SITE 3/8"



Size	Straight tip	Right angled 90° tip
36 Fr	01V362L8	01V362L9

NON-REINFORCED



Size	Straight tip
32 Fr	VC3208