


# Ostial perfusion cannulae

## CORONARY ARTERY OSTIAL CANNULA WITH BALLOON



Balloon size	Tip size	Curved	Straight
4 mm	2,1 mm	OC04A	OC04S
5 mm	2,1 mm	OC05A	OC05S
6 mm	3 mm	OC06A	OC06S
7 mm	3 mm	OC07A	OC07S
8 mm	3,5 mm	OC08A	OC08S


### PRODUCT INFORMATION

Andocor ostial perfusion cannulae are intended for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

**Q** Connection: Female luer lock  
 Length: 30 cm  
 Qty/box: 20 Pcs

**+** Self-inflating silicone balloon  
 Female luer lock connection

### CORONARY ARTERY OSTIAL CANNULA PRESSURE DROP VS FLOW



# Ostial perfusion cannulae

CORONARY ARTERY OSTIAL CANNULA, BASKET TIP WITH FLANGE



Tip size	45° angled	90° angled
10 Fr	OC1045	OC1090
12 Fr	OC1245	OC1290
14 Fr	OC1445	OC1490

## PRODUCT INFORMATION

Andocor ostial perfusion cannulae are intended for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

**Q** Connection: Female luer lock  
 Length: 21 cm  
 Qty/box: 20 Pcs

**+** Soft basket tip  
 Different configurations  
 Stainless steel shaft

## CORONARY ARTERY OSTIAL CANNULA PRESSURE DROP VS FLOW

