

POROUS VENT PLUGS



Connection	QTY/Box	Reference
1/4"	40	VP14
3/8"	40	VP38

VENTED CONNECTOR CAPS



Connection	QTY/Box	Reference
3/8"	40	VCC01

VESSEL CANNULAE



Description	Configuration	Length	Standard
Transparent flexible	Standard	5 cm	VC3
	Closed back check valve	6,5 cm	VC3B
White rigid	Standard	5 cm	VC04

PRODUCT INFORMATION

Andocor vessel cannulae are used to help check for leaks in a harvested vein that will be used in a graft with sterile saline solution. Additionally, the vessel cannula can be used for antegrade cardioplegia administration.

Q Connection: Female luer connector
Qty/box: 40 Pcs

QUICK PRIME LINES



Length	QTY/Box	Reference
100 cm	20	QPL100

PRODUCT INFORMATION

Andocor quick prime lines facilitate the transfer of fluids during the priming of the extracorporeal circuit.

