

Paediatric arterial cannulae

REINFORCED WITH 90° METAL TIP



PRODUCT INFORMATION

BODY LENGTH: 23 CM

Andocor paediatric arterial cannulae are intended for perfusion of the ascending aorta during cardiopulmonary bypass surgery up to six hours or less.

Q Connection: 3/16" - 1/4"
Qty/box: 20 Pcs

+ Kink resistant wire reinforced PVC tubing
Steel tip with thin wall for high flow rates

Description	Size	With vent plug
90° metal tip	4,5 Fr	AP04519V
	6 Fr	AP06019V
	7,5 Fr	AP07519V
	9 Fr	AP09019V
	10,5 Fr	AP10519V
	12 Fr	AP12019V

METAL TIP ARTERIAL CANNULA PRESSURE DROP VS FLOW

