

# Paediatric venous cannulae

REINFORCED WITH 90° METAL TIP



## PRODUCT INFORMATION

BODY LENGTH: 35 CM

Andocor paediatric venous cannulae are intended for drainage from the superior and inferior vena cava during cardiopulmonary bypass surgery up to six hours or less.

**Q** Connection: 10-12 Fr: 3/16" - 1/4" | Other: 1/4"  
Qty/box: 20 Pcs

**+** Kink resistant wire reinforced PVC tubing  
Steel tip with thin wall for high flow rates  
Additional side holes for improved drainage

Description	Size	Reference
90° metal tip	10 Fr	01V101L7
	12 Fr	01V121L7
	14 Fr	01V141L7
	16 Fr	01V161L7
	18 Fr	01V181L7
	20 Fr	01V201L7

## METAL TIP VENOUS CANNULA PRESSURE DROP VS FLOW

