

Pericardial sumps



Type A



Type B



Type C

Description	Size	Tubing	Tip length	Tubing length	Reference
Type A Metal tip	18 Fr	Flexible	2 cm	30 cm	PS40M
	18 Fr	Hard	2 cm	30 cm	PS40MH
	18 Fr	Flexible	2 cm	79 cm	PS40ML
	18 Fr	Hard	2 cm	79 cm	PS40MLH
Type B Perforated pool tip	23 Fr	Flexible	3,5 cm	30 cm	PS40M1
Type C Perforated pool tip	15 Fr	Flexible	2 cm	30 cm	PS40M2
	15 Fr	Hard	2 cm	79 cm	PS40MLH2

PRODUCT INFORMATION

Andocor suckers are used for the removal of excess fluid from the surgical field.



Connection: 1/4"
Qty/box: 20 Pcs



Smooth design
Weighted tip
Stainless steel tip

WIREWOUND TIP

Description	Size	Tubing	Tip length	Tubing length	Reference
Wirewound tip	16 Fr	Flexible	4 cm	30 cm	PS40

