

# Retrograde cardioplegia cannulae



Size	Balloon diameter	Reference
14 Fr	18 mm	RCC01

## PRODUCT INFORMATION

Andocor retrograde cardioplegia cannulae are intended to be connected to the cardioplegia line in order to infuse cardioplegic solution and blood into the patient's coronary sinus during open heart procedures.

**Q** Length: 31 cm  
Qty/box: 20 Pcs

**+** Self-inflating  
Silicone balloon  
Malleable stylet

## RETROGRADE CARDIOPLEGIA CANNULA PRESSURE DROP VS FLOW

