

# Silicone vent catheters



Size	Tip perforation size	Length	Reference
16 Fr	6,5 cm	38 cm	SV16
20 Fr	7,2 cm	39 cm	SV20

## PRODUCT INFORMATION

Andocor vent catheters are medical devices to decompress or vent the heart during cardiopulmonary bypass cannulation. Vent catheters are used as a suction catheter to assist in the complete evacuation of blood.



Connection: 1/4"  
Qty/box: 20 Pcs



Guidewire stylet  
Soft and flexible silicone tubing  
Soft tip