


Thoracic drains

HARD TUBING



Description	Size	Without connector	With connector
Right angled	20 Fr	TC20B	TC20BC
	24 Fr	TC24B	TC24BC
	28 Fr	TC28B	TC28BC
	32 Fr	TC32B	TC32BC
Straight	20 Fr	TC20S	TC20SC
	24 Fr	TC24S	TC24SC
	28 Fr	TC28S	TC28SC
	32 Fr	TC32S	TC32SC

PRODUCT INFORMATION


Andocor thoracic drainage catheters are intended for the removal of air or fluid from the intrathoracic space. The smooth bevelled connection site facilitates the insertion. They are constructed of a clear and kink-resistant PVC tubing. Can be provided with a matching size connector for connection to the drainage system.

Q Length: 58 cm
 Connection: For 20-24 Fr 1/4", for 28-32 Fr 3/8"
 Qty/box: 20 Pcs

+ Kink resistant PVC tubing
 Side holes and open rounded tip on distal end
 Depth markings and radiopaque stripe

Thoracic drains

SOFT TUBING



Description	Size	With connector
Right angled	20 Fr	TCS20B
	24 Fr	TCS24B
	28 Fr	TCS28B
	32 Fr	TCS32B
Straight	20 Fr	TCS20S
	24 Fr	TCS24S
	28 Fr	TCS28S
	32 Fr	TCS32S

PRODUCT INFORMATION

Andocor thoracic drainage catheters are intended for the removal of air or fluid from the intrathoracic space. The tapered tongue at the proximal end facilitates the insertion. They are constructed of a frosted and kink-resistant PVC tubing. Provided with matching size connector for connection to the drainage system.

Q Length: 45 cm
Qty/box: 20 Pcs

+ Kink resistant PVC tubing
Side holes and open rounded tip on distal end
Depth markings and radiopaque stripe