


# Tourniquet sets



## CANNULATION TOURNIQUET SET

| Description                                                   | Snare | Red tube | Blue tube | Clear tube | Reference |
|---------------------------------------------------------------|-------|----------|-----------|------------|-----------|
| Tube size: 11 Fr<br>Tube length: 17 cm<br>Snare length: 25 cm | 1     | 1        | 1         | 3          | CT01      |
|                                                               | 1     | 1        | 2         | 3          | CT02      |
|                                                               | 2     | 1        | 2         | 5          | CT03      |
|                                                               | 1     | 1        | 1         | -          | CT04      |
|                                                               | 1     | 2        | 1         | -          | CT05      |
|                                                               | 2     | 2        | 2         | 4          | CT06      |
|                                                               | 2     | 4        | 2         | 2          | CT07      |
|                                                               | 1     | 5        | -         | -          | CT08      |
|                                                               | 1     | -        | -         | 3          | CT09      |

## PRODUCT INFORMATION

Andocor cannulation tourniquet sets are used to secure arterial and venous cannulae during cardiopulmonary bypass.

-  Tubing: flexible PVC  
Qty/box: 20 Pcs

|                                      |            |      |
|--------------------------------------|------------|------|
|                                      | Clear Tube |      |
| Tube size: 21 Fr, Tube length: 20 cm | 2          | CT10 |