


Flex line two stage venous cannulae

PROXIMAL AND DISTAL REINFORCED WITH LIGHTHOUSE TIP *



ROUND BODY



FLAT BODY



ROUND BODY					FLAT BODY		
Size	Length	Without connector	With connector	With luer connector	Without connector	With connector	With luer connector
32/40 Fr	36 cm	01V32L40	01V32L401	01V32L402	01V32L40FB	01V32L401FB	01V32L402FB
34/40 Fr	36 cm	01V34L40	01V34L401	01V34L402	01V34L40FB	01V34L401FB	01V34L402FB
34/46 Fr	40 cm	01V34L46	01V34L461	01V34L462	01V34L46FB	01V34L461FB	01V34L462FB
36/46 Fr	40 cm	01V36L46	01V36L461	01V36L462	01V36L46FB	01V36L461FB	01V36L462FB
36/51 Fr	40 cm	01V36L51	01V36L511	01V36L512	01V36L51FB	01V36L511FB	01V36L512FB

PRODUCT INFORMATION


Andocor two stage venous cannulae are intended for use in venous drainage via the right atrium and inferior vena cava simultaneously during cardiopulmonary bypass surgery up to six hours or less.

- 

Connection site: 1/2"
Qty/box: 20 Pcs
- 

Single piece construction
Open and flexible lighthouse tip
Easy positioning
Insertion depth marker

VENOUS CANNULA PRESSURE DROP VS FLOW



Special configurations

CONNECTION SITE: 3/8" *



Size	Length	With connector
32/40 Fr	36 cm	01V32L403
34/46 Fr	40 cm	01V34L463

* This product is currently not CE-marked and is not available for commercial distribution within the European Union, unless a derogation applies. Market placement will be possible upon successful completion of the conformity assessment in accordance with Regulation (EU) 2017/745.

For further information, please contact regulatory@andocor.com.