


Ostial perfusion cannulae

WITH BALLOON



Balloon size	Tip size	Curved	Straight
6 mm	3 mm	OC06A	OC06S
7 mm	3 mm	OC07A	OC07S


PRODUCT INFORMATION

Andocor ostial perfusion cannulae are intended for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

Q Connection: Female luer lock
Length: 30 cm
Qty/box: 20 Pcs


+ Self-inflating silicone balloon
Female luer lock connection

CORONARY ARTERY OSTIAL CANNULA
PRESSURE DROP VS FLOW



Ostial perfusion cannulae

BASKET TIP WITH FLANGE



Tip size	45° angled	90° angled
10 Fr	OC1045	OC1090
12 Fr	OC1245	OC1290
14 Fr	OC1445	OC1490

PRODUCT INFORMATION

Andocor ostial perfusion cannulae are intended for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

Q Connection: Female luer lock
Length: 21 cm
Qty/box: 20 Pcs

+ Soft basket tip
Different configurations
Stainless steel shaft

CORONARY ARTERY OSTIAL CANNULA
PRESSURE DROP VS FLOW

