


# SlimFlex arterial cannulae

REINFORCED



BODY LENGTH: 22 CM

| Description         |  | Size  | With blunt introducer | With introducer |
|---------------------|--|-------|-----------------------|-----------------|
| With connector      |  | 18 Fr | A18101P               | A18101PW        |
|                     |  | 20 Fr | A20101P               | A20101PW        |
|                     |  | 22 Fr | A22101P               | A22101PW        |
| With luer connector |  | 18 Fr | A18102P               | A18102PW        |
|                     |  | 20 Fr | A20102P               | A20102PW        |
|                     |  | 22 Fr | A22102P               | A22102PW        |

## PRODUCT INFORMATION \*


This cannula is intended for arterial cannulation and perfusion during cardiopulmonary bypass surgery, up to six hours or less.

The cannula consists of soft PVC tubing with encapsulated stainless steel wire reinforcement. Each cannula is supplied with an introducer, which is designed to accommodate a 0.038" (0.97mm) diameter guidewire.

 Connection site: 3/8"  
Qty/box: 20 Pcs

 Soft PVC reinforced tubing

## SLIMFLEX ARTERIAL CANNULAE



\* This product is currently not CE-marked and is not available for commercial distribution within the European Union, unless a derogation applies. Market placement will be possible upon successful completion of the conformity assessment in accordance with Regulation (EU) 2017/745.

For further information, please contact [regulatory@andocor.com](mailto:regulatory@andocor.com).