



VESSEL CANNULAE

Description	Configuration	Length	Standard
Transparent flexible	Standard	5 cm	VC3
	Closed back check valve	6,5 cm	VC3B
White rigid	Standard	5 cm	VC04

PRODUCT INFORMATION

Andocor vessel cannulae are used to help check for leaks in a harvested vein that will be used in a graft with sterile saline solution. Additionally, the vessel cannula can be used for antegrade cardioplegia administration.

-  Connection: Female luer connector
Qty/box: 40 Pcs