



YANKAUER SUCTION TUBES *

Tubing	Length	Open tip 22 Fr	Bulb tip 30 Fr
Vented	25 cm	ST04	ST02
Non-vented	25 cm	ST03	ST01

PRODUCT INFORMATION

Andocor yankauer suction tubes are intended to aspirate body waste fluid during surgery.

Q Connection: 1/4"
Qty/box: 20 Pcs

* This product is currently not CE-marked and is not available for commercial distribution within the European Union, unless a derogation applies. Market placement will be possible upon successful completion of the conformity assessment in accordance with Regulation (EU) 2017/745.

For further information, please contact regulatory@andocor.com.